Bioinformatics file formats

Jessen V. Bredeson UC Berkeley

Goals and outline

- Understand importance of standardized file formats
- Introduce you to commonly-used formats in bioinformatics
- Resources for manipulating or parsing them yourself

Why are (standardized) file formats important?


Syntactic and semantic interoperability

"The capability to communicate, execute programs, or transfer data among various functional units in a manner that requires the user to have little or no knowledge of the unique characteristics of those units"¹

"the capability of different programs to exchange data via a common set of exchange formats, to read and write the same file formats, and to use the same protocols.... the lack of interoperability can be a consequence of a lack of attention to standardization during the design of a program"²


Data sharing and collaboration

File standards provide:
a common language for data sharing,
promote collaboration,
ensure data reusability,
reduce semantic and syntactic errors



- 1. ISO/IEC 2382-01 Information Technology Vocabulary, Fundamental Terms
- 2. Gordon and Hernandez, The Official Guide to the SSCP Book


Variant-calling workflow



We have a specimen of interest...



We extract DNA...



We have data, now what?

Index the genome

```
$ bwa index Trex_genome.fasta # creates the following files:
$ ls -1 Trex_genome.fasta*
Trex_genome.fasta
Trex_genome.fasta.bwt
Trex_genome.fasta.pac
Trex_genome.fasta.ann
Trex_genome.fasta.amb
Trex_genome.fasta.sa
$ samtools faidx Trex_genome.fasta # creates Trex_genome.fasta.fai
$ head -3 Trex_genome.fasta.fai
Chr1 217471166 141
                           100 101
Chr2 181034961 219646160 100 101
Chr3 153873357 402491612 100 101
$ samtools dict Trex_genome.fasta >Trex_genome.dict
$ head -3 Trex_genome.dict
@HD VN:1.0 SO:unsorted
@SQ SN:Chr1 LN:217471166 M5:56d95ce6647ea9087b857b1efa6d00dd
@SQ SN:Chr2 LN:181034961 M5:20852c561ea38c67aa67e6d655cfebf2
```

FASTA/Pearson

https://en.wikipedia.org/wiki/FASTA_format

>U31202.1 Human noggin (NOGGIN) gene, complete cds GAGCTCCGGCGGGTCAGCCGGACTGTCGGCTTCCCGGGGCATCTGGGTCCGGCGGGGCACAGCCCTGGGC GCTGCCGAAGCCGCCGCCGCCGCCTCCGCGGCGAGTACAGGCGGCTTCCCCGGAGCCTGTGCAGCTCCA GAGAGAGTCAGTGGTTTCCATGGTGATGGAGCTGAAAGTGCAGGAAATTTAAAGGCTTGGACCCTGCGAG ACAGACAAACCGGTGCCAACGTGCGCGGACGCCGCCGCCGCCGCCGCTGGAGTCCGCCGGGCAGAGC AGCGGCCGGNCGAAGAGCAGCGAGAGGAGGAGGGGAGAGCGGCTCGTCCACGCGCCCTGCGCCGCCGCCG GCCCGGGAAGGCAGCGAGGAGCCGGCGCCTCCCGCGCCCCGCGGTCGCCCTGGAGTAATTTCGGATGCCC AGCCGCGGCCGCCTTCCCCAGTAGACCCGGGAGAGGAGTTGCGGCCAACTTGTGTGCCTTTCTTCCGCCC CGGTGGGAGCCGCCGCGAAGGGCTCTCCCGGCGGCTCATGCTGCCGGCCCTGCGCCTGCCCAGCC GACGCGGGACGAAGCAGCACCCGGGCGCGCGCCAGAGGCATGGAGCGCTGCCCCAGCCTAGGGGTCAC AGGAAAAGGATCTGAACGAGACGCTGCTGCGCTCGCTCGGGGGGCCACTACGACCCAGGCTTCATGGC AGGGCTTGGCCCAGGGCAAGAAGCAGCGCCTAAGCAAGAAGCTGCGGAGGAAGTTACAGATGTGGCTGTG GTCGCAGACATTCTGCCCCGTGCTGTACGCGTGGAACGACCTGGGCAGCCGCTTTTGGCCGCGCTACGTG AAGGTGGGCAGCTGCTTCAGTAAGCGCTCGTGCTCCGTGCCCGAGGGCATGGTGTGCAAGCCGTCCAAGT $\mathsf{CCGTGCACCTCACGGTGCTGCGGTGGCGCTGTCAGCGGCGGGGGGCCAGCGCTGCGGCTGGATTCCCAT$ ACTTGATCCTCGAGCTC

FASTA/Pearson

FASTA Defline: Sequence ID + Description on same line, sequence string on the next

>U31202.1 Human noggin (NOGGIN) gene, complete cds

"greater than"
Start of record

Sequence ID

Required;
Any printable
non-whitespace
characters:

[!-~]

Whitespace only required if description present

GASCTCCGGCGGGTCAGCCGGACTGTCGGCTTCCCGGGGGCATCTGGGTCCGGCGGGGCACAGCCCTGGGC **&**CTGCCGAAGCCGCCGCCGCCCCCCGCGGCGAGTACAGGCGGCTTCCCCCGGAGCCTGTGCAGCTCCA GAGAGAGTCAGTGGTTTCCATGGTGATGGAGCTGAAAGTGCAGGAAATTTAAAGGCTTGGACCCTGCGAG ACAGACAAACCGGTGCCAACGTGCGCGGACGCCGCCGCCGCCGCCGCCGCTGGAGTCCGCCGGCAGAGC AGCGGCCGGNCGAAGAGCAGCGAGAGGAGGAGGGGGAGAGCGGCTCGTCCACGCGCCCTGCGCCGCCGCCG GCCCGGGAAGGCAGCGAGGAGCCGCCCTCCCGCGCCCCGCGGTCGCCCTGGAGTAATTTCGGATGCCC AGCCGCGGCCGCCTTCCCCAGTAGACCCGGGAGAGGAGTTGCGGCCAACTTGTGTGCCTTTCTTCCGCCC CGGTGGGAGCCGCCGCGAAGGGCTCTCCCGGCGGCTCATGCTGCCGGCCCTGCGCCTGCCCAGCC GACGCGGGACGAAGCAGCCCCGGGCGCGCCAGAGGCATGGAGCGCTGCCCCAGCCTAGGGGTCAC CGCCCGGCACCCAGCGACAACCTGCCCCTGGTGGACCTCATCGAACACCCAGACCCTATCTTTGACCCCA AGGAAAAGGATCTGAACGAGACGCTGCTGCGCTCGCTCGGGGGGCCACTACGACCCAGGCTTCATGGC AGGGCTTGGCCCAGGGCAAGAAGCAGCGCCTAAGCAAGAAGCTGCGGAGGAAGTTACAGATGTGGCTGTG GTCGCAGACATTCTGCCCCGTGCTGTACGCGTGGAACGACCTGGGCAGCCGCTTTTGGCCGCGCTACGTG AAGGTGGGCAGCTGCTTCAGTAAGCGCTCGTGCTCCGTGCCCGAGGGCATGGTGTGCAAGCCGTCCAAGT ${\sf CCGTGCACCTCACGGTGCTGCGGTGGCGCTGTCAGCGGCGGGGGGCCAGCGCTGCGGCTGGATTCCCAT}$ CCAGTACCCCATCATTTCCGAGTGCAAGTGCTCGTGCTAGAACTCGGGGGCCCCCTGCCCGCACCCGGAC

Description/Comment
CTGGGC
GCTCCA
Optional;
AACCCC
Free-form text

FASTA Body/ Sequence string

> Nucleotide, amino acid, IUPAC codes, alignment characters [-*]

May/may not be wrapped flush

>lcllBC064885.2_cds_AAH64885.1_1 [gene=mtpn] [protein=myotrophin] [protein_id=AAH64885.1 ATGGGTGACAAGGAGTTCGTGGGCCATCAAGAACGGAGACCTGGATGCAGTGAAAGAATTCGTACTTG

CCCCTACTATCTGCCTGCTACGAGGGCCATCGCAAATGTGCTGCTGCTGAAAGATGGCATCACC

FASTA files are best suffixed with ".fasta" or ".fa"; some tools require this.

CCATTAA

ACTTGATCCTCGAGCTC

https://en.wikipedia.org/wiki/FASTQ_format

FASTQ Sequence Header: Sequence ID + Description on same line, sequence string on the next

, Whitespace only required if description present

"At" symbol

Start of sequence portion of record

Sequence ID

Required;
Any printable
non-whitespace
characters
[!-~]

Description/Comment optional

GGATCTATGGCCATGTAGGGACCATCTGAAGGCAGATCAAAATTTCGCTGAGCAAATTTAGGGTCCGGGTTTGTT

ATAAAAAAAAATTAATTATTTTAAAACTAATTTTTAAATTAATTGGTTTTTTGTGGAATGGTATT

+

AAFFFJFJAJJJFFA<JF-7FJF<JJJJ--<<FJ-J<A7---<-FFJFJFJJAJ-F<F<<-F-7---7-<-<FF @SRR10178655.3 0:N:0:

CATTATATACGTCGCCACTCTTAATTTCCTTTTCCATAAGAGCGTATAATCTTGTAATACAATGTCTTCTCCAA

+

AAAFFJJJJAFJJJJAFJF</ri>
AAAFFJJJJF
AAAFFJJJJF
AAAFFJJJJF
AAAFFJJJJAFJJJF
AAAFFJJJJF
AAAFFJJJJAFJJJF
AAAFFJJJJF
AAAFFJJJJJF
AAAFFJJJJF
AAAFFJJJJF
AAAFFJJJJJF
AAAFFJJJJF
AAAFFJJJJF
AAAFFJJJF
AAAFFJJJF
AAAFFJJJF
AAAFFJJJF
AAAFFJJJF
AAAFFJJJF
AAAFFJJF
<pre

AAGTTATTCTGCCTCTAATGCGATAACTGTAATCTTTAATTGTGTAATTTCTTTTTCACAATCTGAGCCACGCCA

AAAAAJF<-A<-7FJFJJJJFJFJJJJJ<FJ--7<FF-7-<--7-A<7FFJAJFFJJJAJ7FF-F7FA-7<-A-7-

+

AAFFFJFJJJJJJ<FFJJJAFJJ<FFAJFJ--FJJJJ-FFF-<FFJFA-FJJ-AJ-<<-FFFAFJJJJJJAJ-7--

FASTQ Sequence

Nucleotide, amino acid, IUPAC codes

Typically *not* wrapped flush (discouraged)

FASTQ files are best suffixed with ".fastq" or ".fq", some tools require this.

FASTQ Sequence Header: Paired-end or mate-pair reads

CATTTTTCCAAACATACC

<FFFJFFNJFJJJJ-F-<FJ7JJFJJF<F-7A-7FJ-<FJJ<<FJFFJJFJ<

Read 2

Read 1

FASTQ Qualities Header: (optional) Same as Sequence Header, or absent completely

"Plus" symbol_ Start of qualities portion of record

Qualities ID
Optional;
If present, typically same as Sequence
ID; Must follow same rules

AAFFFJFJJJJJJJ<FFJJAFJJ<FFAJFJ--FJJJJ-FFF-<FFJFA-FJJ-AJ-<<-FFFAFJJJJJJAJ-7---

FASTQ Qualities ASCII+offset encoded "Phred"

scores.

Must be same length as sequence.

Typically *not* wrapped flush (discouraged)

 $Phred = -10 \bullet \log_{10}(P)$

P = fractional probability that the base call is wrong

| ▼
Dec | H | Oct | Chai | r I | ↓
Dec | Нх | Oct | Html | Chr | ↓
Dec | Нх | Oct | Html | Chr | l Dec | : Hx | Oct | Html Ch | hr |
|-----------------|---|-----|------|--------------------------|-----------------|----|-----|---|-------|-----------------|----|-----|-------|-----|-------|------|-----|---------|----|
| 0 | | | | (null) | | | | | Space | | | | @ | | | | | ` | N. |
| 1 | | | | (start of heading) | | | | a#33; | _ | | | | a#65; | | | | | a#97; | a |
| 2 | | | | (start of text) | | | | a#34; | | 66 | | | a#66; | | | | | | b |
| 3 | | | | (end of text) | | | | a#35; | | | | | a#67; | | I | | | | c |
| 4 | | | | (end of transmission) | | | | a#36; | | 68 | | | a#68; | | | | | d | d |
| 5 | | | | (enquiry) | | | | a#37; | | 69 | 45 | 105 | a#69; | E | | | | e | |
| 6 | | | | (acknowledge) | 38 | 26 | 046 | @#38; | 6 | 70 | 46 | 106 | a#70; | F | 102 | 66 | 146 | a#102; | f |
| 7 | 7 | 007 | BEL | (bell) | 39 | 27 | 047 | @#39; | 1 | 71 | 47 | 107 | G | G | 103 | 67 | 147 | @#103; | g |
| 8 | 8 | 010 | BS | (backspace) | 40 | 28 | 050 | &# 4 0; | (| 72 | 48 | 110 | H | H | 104 | 68 | 150 | a#104; | h |
| 9 | 9 | 011 | TAB | (horizontal tab) | 41 | 29 | 051 |) |) | 73 | 49 | 111 | 6#73; | I | 105 | 69 | 151 | i | i |
| 10 | Α | 012 | LF | (NL line feed, new line) | 42 | 2A | 052 | &#42;</td><td>*</td><td>74</td><td>4A</td><td>112</td><td>a#74;</td><td>J</td><td>106</td><td>6A</td><td>152</td><td>j</td><td>j</td></tr><tr><td>11</td><td>В</td><td>013</td><td>VT</td><td>(vertical tab)</td><td>43</td><td>2B</td><td>053</td><td>&#43;</td><td>+</td><td>75</td><td>4B</td><td>113</td><td>a#75;</td><td>K</td><td>107</td><td>6B</td><td>153</td><td>k</td><td>k</td></tr><tr><td>12</td><td>С</td><td>014</td><td>FF</td><td>(NP form feed, new page)</td><td>44</td><td>20</td><td>054</td><td>a#44;</td><td>,</td><td>76</td><td>4C</td><td>114</td><td>a#76;</td><td>L</td><td>108</td><td>6C</td><td>154</td><td>4#108;</td><td>1</td></tr><tr><td>13</td><td>D</td><td>015</td><td>CR</td><td>(carriage return)</td><td>45</td><td>2D</td><td>055</td><td>a#45;</td><td>F 1.</td><td>77</td><td>4D</td><td>115</td><td>@#77;</td><td>M</td><td></td><td></td><td></td><td>m</td><td></td></tr><tr><td>14</td><td>E</td><td>016</td><td>so</td><td>(shift out)</td><td>46</td><td>2E</td><td>056</td><td>&#46;</td><td></td><td>78</td><td>4E</td><td>116</td><td>۵#78;</td><td>N</td><td>110</td><td>6E</td><td>156</td><td>@#110;</td><td>n</td></tr><tr><td>15</td><td>F</td><td>017</td><td>SI</td><td>(shift in)</td><td></td><td></td><td></td><td>6#47;</td><td></td><td>79</td><td>4F</td><td>117</td><td>%#79;</td><td>0</td><td></td><td></td><td></td><td>o</td><td></td></tr><tr><td>16</td><td>10</td><td>020</td><td>DLE</td><td>(data link escape)</td><td>48</td><td>30</td><td>060</td><td>@#48;</td><td>0</td><td>80</td><td></td><td></td><td>P</td><td></td><td></td><td></td><td></td><td>p</td><td></td></tr><tr><td></td><td></td><td></td><td></td><td>(device control 1)</td><td></td><td></td><td></td><td>a#49;</td><td></td><td></td><td></td><td></td><td>Q</td><td></td><td></td><td></td><td></td><td>q</td><td></td></tr><tr><td>18</td><td>12</td><td>022</td><td>DC2</td><td>(device control 2)</td><td>50</td><td>32</td><td>062</td><td>2</td><td>2</td><td></td><td></td><td></td><td>4#82;</td><td></td><td>114</td><td>72</td><td>162</td><td>r</td><td>r</td></tr><tr><td>19</td><td>13</td><td>023</td><td>DC3</td><td>(device control 3)</td><td></td><td></td><td></td><td>3</td><td></td><td></td><td></td><td></td><td>S</td><td></td><td></td><td></td><td></td><td>s</td><td></td></tr><tr><td>20</td><td>14</td><td>024</td><td>DC4</td><td>(device control 4)</td><td></td><td></td><td></td><td>4</td><td></td><td></td><td></td><td></td><td>۵#84;</td><td></td><td></td><td></td><td></td><td>t</td><td></td></tr><tr><td>21</td><td>15</td><td>025</td><td>NAK</td><td>(negative acknowledge)</td><td></td><td></td><td></td><td>5</td><td></td><td></td><td></td><td></td><td><u>4</u>#85;</td><td></td><td>I — — ·</td><td></td><td></td><td>u</td><td></td></tr><tr><td></td><td></td><td></td><td></td><td>(synchronous idle)</td><td></td><td></td><td></td><td>a#54;</td><td></td><td>86</td><td></td><td></td><td>4#86;</td><td></td><td>I — — -</td><td></td><td></td><td>v</td><td></td></tr><tr><td></td><td></td><td></td><td></td><td>(end of trans. block)</td><td></td><td></td><td></td><td><u>@</u>#55;</td><td></td><td>87</td><td></td><td></td><td><u>6#87;</u></td><td></td><td></td><td></td><td></td><td>w</td><td></td></tr><tr><td></td><td></td><td></td><td></td><td>(cancel)</td><td></td><td></td><td></td><td>a#56;</td><td></td><td>88</td><td></td><td></td><td>X</td><td></td><td></td><td></td><td></td><td>x</td><td></td></tr><tr><td></td><td></td><td>031</td><td></td><td>(end of medium)</td><td>57</td><td></td><td></td><td>a#57;</td><td></td><td>89</td><td></td><td></td><td>6#89;</td><td></td><td>ı</td><td></td><td></td><td>y</td><td></td></tr><tr><td></td><td></td><td>032</td><td></td><td>(substitute)</td><td>58</td><td></td><td></td><td>6#58;</td><td></td><td></td><td></td><td></td><td>6#90;</td><td></td><td>1</td><td></td><td></td><td>z</td><td></td></tr><tr><td></td><td></td><td>033</td><td></td><td>(escape)</td><td></td><td></td><td></td><td>6#59;</td><td></td><td>91</td><td></td><td></td><td>6#91;</td><td>-</td><td>123</td><td></td><td></td><td>{</td><td></td></tr><tr><td></td><td></td><td>034</td><td></td><td>(file separator)</td><td></td><td></td><td></td><td>4#60;</td><td></td><td></td><td></td><td></td><td>6#92;</td><td></td><td>I — — -</td><td></td><td></td><td> </td><td></td></tr><tr><td></td><td></td><td>035</td><td></td><td>(group separator)</td><td></td><td></td><td></td><td>=</td><td></td><td></td><td></td><td></td><td>6#93;</td><td>-</td><td>I — — -</td><td>. –</td><td></td><td>}</td><td></td></tr><tr><td></td><td></td><td>036</td><td></td><td>(record separator)</td><td>62</td><td></td><td></td><td>></td><td></td><td></td><td></td><td></td><td>a#94;</td><td>^</td><td></td><td></td><td></td><td>~</td><td></td></tr><tr><td>31</td><td>1F</td><td>037</td><td>US</td><td>(unit separator)</td><td>63</td><td>ЗF</td><td>077</td><td>?</td><td>2</td><td>95</td><td>5F</td><td>137</td><td>_</td><td>_</td><td> 127</td><td>7F</td><td>177</td><td></td><td>DEL</td></tr></tbody></table> | | | | | | | | | | | |

| P | Phred |
|--------------------|-------|
| 1×10 ⁰ | 0 |
| 1×10 ⁻¹ | 10 |
| 1×10 ⁻² | 20 |
| 1×10 ⁻³ | 30 |
| 1×10 ⁻⁴ | 40 |
| 1×10 ⁻⁵ | 50 |
| 1×10 ⁻⁶ | 60 |

Source: www.LookupTables.com


 $Phred = -10 \bullet \log_{10}(P)$

P = fractional probability that the base call is wrong

```
!"#$%&'()*+,-./0123456789:;<=>?@ABCDEFGHIJKLMNOPQRSTUVWXYZ[\]^_`abcdefghijklmnopqrstuvwxyz{|}~
33
            59
               64
                   73
                                            126
                                 104
0.....9.......40
                S - Sanger
        Phred+33, raw reads typically (0, 40)
        Solexa+64, raw reads typically (-5, 40)
X - Solexa
I - Illumina 1.3+ Phred+64, raw reads typically (0, 40)
J - Illumina 1.5+ Phred+64, raw reads typically (3, 41)
  with 0=unused, 1=unused, 2=Read Segment Quality Control Indicator (bold)
  (Note: See discussion above).
L - Illumina 1.8+ Phred+33, raw reads typically (0, 41)
```

QA/QC'ing Illumina Reads

```
$ fastqc --threads 2 --extract SRR10178655_1.fastq.gz SRR10178655_2.fastq.gz
# View the FastQC results in Safari (Mac only):
$ open -a Safari.app SRR10178655_1_fastqc/fastqc_report.html
```



Trim and align Reads

- \$ java -Xmx500m -jar ./Trimmomatic-0.39/trimmomatic-0.39.jar PE -phred33 \
 -summary SRR10178655.summary SRR10178655_1.fastq.gz SRR178655_2.fastq.gz \
 SRR10178655_1_passed.fastq.gz SRR10178655_1_failed.fastq.gz \
 SRR10178655_1_passed.fastq.gz SRR10178655_2_failed.fastq.gz MINLEN:100 \
 ILLUMINACLIP:./Trimmomatic-0.39/adapters/NexteraPE-PE.fa:2:30:10:2:keepBothReads
- \$ bwa mem -R '@RG\tID:SRR10178655\tSM:Trex\tLB:HAMMOND01\tPL:ILLUMINA' \
 Trex_genome.fasta SRR10178655_1_passed.fastq.gz SRR10178655_2_passed.fastq.gz | \
 samtools view -b >SRR10178655.bam
- \$ samtools sort -m 1g -o SRR10178655.srt.bam SRR10178655.bam
- \$ samtools index SRR10178655.srt.bam # creates SRR10178655.srt.bam.bai

http://samtools.github.io/hts-specs/SAMv1.pdf http://samtools.github.io/hts-specs/SAMtags.pdf

SAM: Sequence Alignment/Map format

BAM: Binary SAM


CRAM: Reference-Compressed SAM (also binary)

```
@HD VN:1.3 S0:coordinate
@SQ SN:ref|NC_001133| LN:230218
@SQ SN:ref|NC_001134| LN:813184
@SQ SN:ref|NC_001148| LN:948066
@SQ SN:ref|NC_001224| LN:85779
@PG ID:bwa PN:bwa VN:0.7.15-r1140 CL:bwa mem ...
@RG ID:SRR10178655 SM:Trex LB:HAMMOND01 PL:ILLUMINA
SRR10178655.85923 163 ref|NC_001133| 1 30 2S7M = 383 392
    ACATTACTC AAA))-*## NM:i:0 MD:i:7 AS:i:7 RG:Z:SRR10178655
SRR10178655.85923 83 ref|NC_001133| 383 60 9M = 1 -392
    ACCTCACAT 7JFFFFFAA NM:i:0 MD:7:9 AS:i:9 RG:7:SRR10178655
```

SAM Header: Meta information describing file format and data within. Header lines must start with "@" symbol (and read IDs must not). Tab separated. Reference IDs cannot be "*", "0", or "="; they have special meaning.

```
Header format version and sort order
                                                                         Read Group
                                                                       Almost required;
                                                      Program
                    VN:1.3 SO:coordinate
                                                                       ID, sample name,
                                                 processing history
                    SN:ref|NC_001133| LN:230218
              @SO
                                                                       and library names,
                                                 (with commands)
              @SQ
                   SN:ref|NC_001134|
                                      LN:813184
                                                                      sequencing platform
              @SO
                   SN:ref|NC 001148|
                                      LN:948066
 Sequence
                   SN:ref|NC_001224|
              @SQ
                                      LN:85779
 Reference
              @PG
                   ID:bwa PN:bwa
                                   VN:0.7.15-r1140 CL:bwa mem
sequence IDs
              @RG
                   ID: SRR10178655
                                   SM:Trex
                                            LB: HAMMONDØ1
and lengths;
                                      refINC 001133
               SRR10178655.85923
                                163
listed in same
                ACATTACTC AAA))-*##
 order as in
                                       refINC 0011331
               SRR10178655, 85923
                 ACCTCACAT
                                      NM:i:0
                                              MD:7:9
   FASTA
```

SAM Body: Describes mapping and alignment without the reference. Eleven required fields. Tab separated. Undefined values: "0" for numeric field, a "*" for non-numeric.



CIGAR operators

Op Meaning

M : Match

I : Insertion

D : Deletion

= : Sequence match

X : Sequence mismatch

N : Forward-skip query on reference (intron)

H : Query hard clipping

S : Query soft clipping

P : Padded reference

B : Backward-skip query on reference

Example:

Q: ATGACAGGACAGAT-GA^{GG}

R: ATG-CAGGCCAGATTGATA

3M 1I 10M 1D 2S describes same alignment as 3= 1I 4= 1X 5= 1D 2S

Bit Flags

To add or test for flags, use 2ⁿ values with bitwise operations:

```
Add flag(s) Test for flag(s) flags l=2**0 flags & 1024 # correct flags l=2**1 flags l=2**6
```

QA/QC'ing Alignments

```
$ gatk MarkDuplicates --java-options '-Xmx1G' \
    -MAX_FILE_HANDLES 2000 \
    -I SRR10178655.srt.bam \
    -0 SRR10178655.srt.mdup.bam \
    -M SRR10178655.metrics


$ samtools stats -@4 --ref-seq Trex_genome.fasta \
    SRR10178655.srt.mdup.bam >SRR10178655.stats

$ plot-bamstats -s Trex_genome.fasta >Trex_genome.gc

$ plot-bamstats -r Trex_genome.gc -p SRR10178655 SRR10178655.stats

$ open -a Safari.app SRR10178655.html
```

QA/QC'ing Alignments



13,797,114,082

12,862,420,009

1.21%

(93.2%)

total: mapped:

error rate:

\$ samtools tview SRR10178655.srt.mdup.bam Trex_genome.fasta

| • • • | 🁚 jessenbredeson — ssh -Y bredeson@cori.nersc.gov — 128×48 |
|---|--|
| | 211 2221 2231 2241 2251 2261 2271 2281 2291 |
| | GGGAGAGAGGGGCGCACAGACAAGGTAGCCTTGCCGGCTAGCAATCCTCAGCGTACTCTACTTTCTGCTGCCTCTGCATTAGCATAGGGAGAGAGA |
| • · · · · · · · · · · · · · · · · · · · | YM |
| <u></u> | |
| | TC |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | g,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | A.CA |
| | |
| | |
| G | |
| ,,,,,,,,,,,g,,,,,,,,,,,,,,,,,,,,, | |
| | <u>.TTC.C</u> .A |
| | |
| CAA | G ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | • • • • • • • • • • • • • • • • • • • |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| ,,,,,,,,,,,g,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,, |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | |
| | C,,,,,,,,,,,,,,,,,,,,,,,,,, |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | |
| | |
| G | TC.C |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| ,,,,,,,,,,,g,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,c,,,,,,,,,,,,,,,,,,,,,,,,,t, |
| G | |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | C |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | |
| G | |
| | |
| | C |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| CG | ,,,,,,c,,,,g __ ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | c,,,,,,,,,,,t,,,,,t,,,,,,,,,,,,,,, |
| | |
| | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | |
| | C |
| | |

- \$ samtools view -b -f3 -F3852 SRR10178655.srt.mdup.bam > SRR10178655.srt.mdup.proper.bam
 \$ samtools index SRR10178655.srt.mdup.proper.bam
 \$ samtools tview SRR10178655.srt.mdup.proper.bam Trex_genome.fasta

| • • • | 🏦 jessenbred | eson — ssh -Y bredes | on@cori.nersc.gov — 1 | 28×48 | |
|--|---|---|---|---|---|
| 2181 2191 2201 CCTCTACTTTCTACTGCCTCTGCATTA | 2211 2221 | 2231 2241 | 2251 2261 | 2271 2281 | 2291 |
| RR. | | | | | |
| | | | | | |
| ,,,,,,,,,,g,,,, | | | | | |
| | ,,c,,,,,,,,,,,c,,, | | ,,,,,,,,,,,t,,,,c | | ,, ,,,,,,,,,c,g,, |
| ,,,,g,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,g,t,,,a,,,,,,, | ,,,,,,,,,,,t,,,,c | | |
| | | | | • | |
| ,,,,,,,,,,,g,,,,,,,,,,,,,,,,,,,,,,,,,, | | ,,,g,,,,,,,,,, | ,,,,,,,,,,,,g,,,,c | | |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | T | | |
| | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| ,,,,,,,,,,g,,,,,,,,,,,,,,,, | | | | • | |
| | | | | | |
| | | | | :,,,,,t,,,,,,,,,,,,,,, | |
| ,,,,,,,,,,g,,,,,,,,,,, | | | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | | | | | |
| | | | | | A |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | | | G | |
| | | | | | A T |
| | | | | • | |
| | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | | ,, ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | | | | | |
| | | | | | |
| | | | | T. ,,,g,,, | |
| ,,,,,,,,,,g,,,,,,,,,,,,,,, | | | | | |
| | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,c,c,g,, |
| | | | | | |
| | | | | | |
| | | | | | .C. ,,,,g,,,,, |
| | | a | ,,,,,,,,,,,,,a,,,,, | | ,,,, |
| | | | | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | | | |
| ,,,,,,,,,,g,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | ,,,,,,,,,,,t,,,,c | | |
| ,,,,,,,,,,,g,,,,,,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | | | ,,,,,,,,,, |
| | | | | ,,,,g,,,,,,,,,,,,,,,, | |
| C | | | | | |
| CC, | ,,,,g,,,,,,,,,c,,,, | | | .,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| | | | | | |
| | ,,,,,t,,,,,,,,,,,,, | | | | |
| | | | ,,,,,,,,,,,,,t,,,,, | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | | ,,,,,,,,,,,c,,,a,,,,, | |
| ,,,,,,,,,cg,,,,,,,,,,,,, | | | | | |
| | | | | | |

\$ samtools tview SRR10178655.srt.mdup.proper.bam Trex_genome.fasta

| • • • | | 🏦 jes | senbredes | on — ssh -\ | / bredeson | @cori.nersc | .gov — 128 | ×48 | | |
|---|--------------------------------|---------------------------------------|-----------|-------------|------------|---|---|---------------------------------------|--|---|
| 2181 2191 | 2201 | | 2221 | 2231 | 2241 | 2251 | 2261 | 2271 | 2281 | 2291 |
| CCTCTACTTTCTACTGCC | ICIGCATTAGCA | | | | | | | | CIGCCICIGCA | TTAGCATAGGGAGAGAG |
| | | | | | | | | | T | |
| ,,,,,,,,,,g,,, | | | | | | | | | | |
| | ,,,, ,,,,C | ,,,,,,,,,, | ,,C,,,,, | ,,,,,,,,,, | | ,,,,,,,,,, | ,t,,,,c,, | ,,,,,,,,,, | | ,,, , ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | , , , , , g , , , , , , | ,,,,,,,,,,, | ,,C,,,,, | ,g,t,,,a,, | | ,,,,,,,,,, | ,t,,,,c,, | , , , , , , , , , , , , , , , , , , , | ,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | | | | | | | | | | |
| CAA | | • | | | | | | | ,,,,,,,,,,,, | |
| | | | | | | | | | | |
| | | ,,,,,,,,,,,, | | g,,,,, | ,,,,,,,,, | ,,,,,,,,,, | ,,,,,,C,, | , , , , , , , , , , | ,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| ,,,,,,,,,,g,,,,g | , , , , , , , , , , , , , | ,,,,,,,,, | | | | | | | | |
| | ,,,,,,,,,,,, | | | | | • | • | | | |
| | | | | | | + | + | + | | |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | C | | | | | | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | , , , , , , , , , , , , , , , , , , , | ,,,,,,,, | ,,,,,,,,, | | ,,,,,,,,,, | | , , , , , , , , , , , , , , , , , , , | | |
| | | | | | | | C | | | A |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | | | | | • | | . G | | |
| ,,,,,,,,, | | | | | A | | C | | | AT |
| | | | | | | • • • • • • • • • • • | • | | | |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,, | ,,,,,,,,,,, | | ,,,,,,,,,, | | ,,,,,,,,,, | | ,,,,,,,, | ,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | | | | | | | | | | |
| | ,,,,,,,,,,,, | | | | | | | | | |
| | | | | | | | | | | ,,,,,,,,,,g,,,, |
| ,,,,,,,,,,g,,,,g | ,,,,,,,,,,,, | , , , , , , , , , , , | ,,C,,,,, | ,,,,,,,,,, | | ,,,,,,,,,, | ,t,,,,c,, | , | | |
| | ,,,,,,,,,,,, | ,,,,,,,,,,, | | ,,,,,,,,,, | ,,,,,,,,, | ,,,,,,,,,, | | , , , , , , , , , , , | ,,,,,,, | ,,,,, <mark>,</mark> ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,, | ,,,,,,,,,,, | ,,,,,,,,, | ,,,,,,,,, | | ,,,,,,,,,,, | 1111111111 | ,,,,,,,,,, | ,,,,,,,, | ,,,,,,,,,c,c,g,, |
| | | | | | | C. | | | | ,,,,,,,t,,, |
| | | | | | | | C | | | · · · · · · · · · · · · · · · · · · · |
| , , , , , , , , , , , , , , , , , , , | | | | | | | | , , , , , , , , , , | | ,, ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| C G | | | ,,c,,,a, | ,,,,,,,,,, | ,,,,,,,,, | ,,,,,,,,,, | ,a,,,,c,, | ,,,,,,,,,, | ,,,,,,,,,,, |
 |
| ,,,,,,,,,,g,,,,g | • | | | | | | | | | |
| | ,,,,,,,,,,,, | | | | | ,,,,,,,,,,, | | | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,, | | | | | | | ,,,,,,,,,, | | |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | ,,,,,,,,,,,, | | | | | | | | | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | AC | C . | | .GGA | G T C | C | | . C | | |
| C | ,,,, | | | | | ,,,,,,,,,, | ,,,,,,C,, | , , , , , , , , , , , , , , , , , , , | ,,,t,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| | | | | A | | | | | | C |
| | ,,,,,,,,,,,, | | ,,,,,,,,, | ,,,,,,,,, | | ,,,,,,,,,, | | ,,,,,,,,,, | | |
| | ,,,,,,,,,,, | ,,t,,,,,,,, | ,,,,,,,,, | | | ,,,,,,,,,, | | | ,,,,,,,,,,, | ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, |
| 0 | | | C | 11111 | | | , T, , , , , C, , | | | |
| ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,, | | | | | | | | | ,,,a,,,,,, | |
| | | | ,,,,,,,,, | | | | , -, , , , , , , , , | | | |

Call Variants

```
gatk HaplotypeCaller \
    --minimum-mapping-quality 30 \
    --min-base-quality-score 20 \
    --read-validation-stringency SILENT \
    --reference Trex_genome.fasta \
    --input SRR10178655.srt.mdup.proper.bam \
    --output SRR10178655.vcf
```

http://samtools.github.io/hts-specs/VCFv4.3.pdf

VCF: Variant Call Format

Chr1

2192815 .

GG

BCF: Binary VCF ##fileformat=VCFv4.2 ##FILTER=<ID=LowQual, Description="Locus is low quality"> ##FILTER=<ID=PASS,Description="Locus passes all filters"> ##FORMAT=<ID=GQ, Number=1, Type=Integer, Description="The Phred-scaled prob. of the genotype"> ##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype"> ##FORMAT=<ID=AD, Number=R, Type=Integer, Description="Number of observation for each allele"> ##INFO=<ID=DP, Number=1, Type=Integer, Description="Total read depth at the locus"> ##contig=<ID=Chr1,length=217471166> ##contig=<ID=Chr2,length=181034961> ID REF ALT #CHROM POS OUAL FILTER INFO FORMAT Trex 8.826 ./.:0:0,1Chr1 534 Т Α LowQual DP=1 GT:GQ:AD Chr1 1315 G 564.103 PASS DP=51 GT:GQ:AD 110:99:26,25 Chr1 CTC CC 209.026 . GT:GQ:AD 0|1:99:19,12 369655 DP=31 912.199 . DP=36 GT:GO:AD 211:43:0,28,8 Chr1 672396 GTT GT,GGT

DP=64

GT:GQ:AD

0/1:99:46,18

GGTATTTTTAG 253.597 .

VCF Metadata Lines: For humans and computers. Required by most tools to pre-declare how to parse file body correctly. **fileformat Meta**FILTER Meta

Required on first line; explicitly defines soft -Tells tools how to interpret rest of file filters one expects to see ##fileformat=VCFv4.2 in the FILTER column ##FILTER=<ID=LowQual, Description="Locus is low quality"> _ ##FILTER=<ID=PASS, Description="Locus passes all filters"> ##FORMAT=<ID=GQ, Number=1, Type=Integer, Description="The Phred-scaled prob. of the genotype"> ##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype"> ##FORMAT=<ID=AD, Number=R, Type=Integer, Description="Number of observation for each allele"> ##INFO=<ID=DP, Number=1, Type=Integer, Description="Total read depth at the locus"> ##contig=<ID=Chr1,length=217471166> ##contia=<ID=Chr2,length=181034961> #CHROM POS OUAL FILTER INFO **FORMAT** DP=1 Chr1 534 8.826 Low0ual GT:GO:AD DP=51 GT: GOFORMAT9Meta25 Chr1 1315 564.103 PASS 209. INFO. Meta DP=31 Chr1 TExplicitly defines the Chr1 Explicitly defines the = 36 672396 types data to be GT, GGT Chr1 types of Key=Value P=64 21 conting Meta observed in sample Optional, encouraged; data to be observed in column(s) Describes reference sequences **INFO** column

observed in CHROM column

VCF Header Line: Defines columns, including the sample names. Required by most tools to parse file correctly; undefined fields set to "."

```
Locus ID
                                                           Locus-Level Meta
                           if applicable
Chromosome name
                                                              Information
                       e.g., DBsnp ID, etc.
                                           is low quality" Key=Value pair info
    and Position
Sequence IDs should cription Locus
                                       Locus-level iteabout the locus (and all
#be in conting Meta; ber=1, Type=Intege Quality: Score "The samples: at the locus) the genotype ">
#PositionsID=basedber=1, Type=String Phred-scaled Genotype">
##FORMATA<ID=AD, Number=R, Type=Integrob ethat locus is umber of object at locus is umber of object at locus is umber of object at locus.
        ID=DP, Number=1, Type=Integernotseally ivariantal read depth at the locus">
       1a=<ID=Chr1.lenath=217471166>
##cortig=xID=Chr2,length=181034961>
                                         QUAL
                                                             INFO
                                                                    FORMAT
#CHROM
                                                   FILTER
                                                                               Trex
                      REF ALT
                                                                               ./.\0:0,1
                                                                    GT:GO:AD
                                                   Low0ual
                                                            DP=1
Chr1
        534
                                                                    (T:G0:AD
Chr1
         1315
                                                   PASS
                                                            DP=31 /GT:GQ:AD
                                                                               0|1:99:19.12
                                   Locus-level
Chr1 Reference and
                                                           Sample-Level AD
                                                                               2 | Sample Field
                                                                               Contains sample
Chr Alternate Alleles G
                            GGTATTTSoft Filter(s)
                                                         Field Formatting
                             "PASS" = passes filters
  Alleles observed in
                                                           Ordered list of
                                                                                 genotype and
                              "." = no filters applied
  reference sequence
                                                          fields present in
                                                                                associated info
                              anything else = failure
  and samples at the
                                                              samples
                                                                                  at the locus
          locus
```

VCF Loci: Tab-delimited columns. Alleles indexed from 0 (REF) to N (ALT) alleles. Genotypes represented with those indices

```
##fileformat=VCFv4.2
##FILTER=<ID=LowQual, Description="Locus is low quality">
##FILTER=<ID=PASS, Description="Locus passes all filters">
##FORMAT=<ID=GO, Number=1, Type=Integer, Description="The Phred-scaled prob. of the genotype">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
*Substitution\locuser=R,TyComplex locusiption="Number of observation for each allele">
##INFO=<ID=DP, Number=1, Type=InMultirablele; iption="Total read depth at the No=call or hard-
##contig=<ID=Chr1,\engtheletion1and substitution!
                                                                         filtered genotype
##contia=<ID=Chr2,length=181034961
                    REF
#CHROM
        POS
                          ALT
                                        OUAL
                                                 FILTER
                                                          INFO
                                                                 FORMAT
                                                                           Trex
                                        8.826
                                                 LowQual
                                                                            ./.:0:0,1
Chr1
        534
                          Α
                                                          DP=1
                                                                 GT:GQ:AD
                          G
                                                                           110:99:26,25
Chr1
        1315
                                        564.103
                                                 PASS
                                                          DP=51
                                                                 GT:GQ:AD
Chr1
        369655
                     CTC
                                        209.026
                                                          DP=31
                                                                 GT:GQ:AD
                                                                           0|1:99:19,12
                                                                           2|1:43:0,28,8
Chr1
        672396
                          GT, GGT
                                        912.199
                                                          DP=36
                                                                 GT:GQ:AD
                     GTT
        2192815
                                                                 GT:GQ:AD
                                                                           0/1:99:46.18
Chr1
                     GG
                          GGTATTTTTAG
                                       253.597
                                                          DP=64
                                               Phased genotypes
Deletion locus
                                                                             Allele Depth
                                                                             Read count for
                                             Unphased genotype
```

each allele

Insertion locus

Annotation files

BED


https://genome.ucsc.edu/FAQ/FAQformat.html#format1

BED: Browser Extensible Data format

```
track itemRgb="On"
Chr1
     0
             126500
                            0 + 0
                                          126500
                     CLEAR1
                                                  0,0,0
     126500
            128500
                     BREAK1
                            0 +
                                 126500
                                          128500
                                                  213,221,213
Chr1
Chr1
     128500
             278000
                     CLEAR2
                            0
                               +
                                 128500
                                          278000
                                                  0,0,0
                     BREAK2
     278000
             280000
                            0 + 278000
                                          280000
                                                 213,221,213
Chr1
Chr1
     280000
             362500
                     CLEAR3
                            0 + 280000
                                          362500
                                                  0,0,0
Chr1
     362500
             366000
                     BREAK3
                            0
                               +
                                 362500
                                          366000
                                                 213,221,213
     366000
Chr1
             427500
                     CLEAR4
                                  366000
                                          427500
                                                  0,0,0
                     BREAK4
     427500
             429500
                               + 427500
                                          429500
                                                 213,221,213
Chr1
                            0
Chr1
     429500
             599500
                     CLEAR5
                               + 429500
                                          599500
                                                  0,0,0
     599500
Chr1
             605500
                     BREAK5
                                  599500
                                          605500
                                                  213,221,213
                            0 +
```

BED

BED: Columns tab-delimited. First three required, all others optional (first 6 typical).



GFF3

https://github.com/The-Sequence-Ontology/Specifications/blob/master/gff3.md

GFF: Generic Feature Format

```
##aff-version 3
##species http://www.ncbi.nlm.nih.gov/Taxonomy/Browser/wwwtax.cgi?id=436495
##genome-build RexBase Trex1
##sequence-region Chr1 1 217471166
# Note Trex_genome.fasta, complete genome
Chr1 Gnomon gene
                                78350
                                                 ID=gene32251; Name=L0C101732307
                          43895
Chr1 Gnomon
             mRNA
                          43895 78350
                                                ID=rna61088; Name=XM_012954515.1; Parent=gene32251
             CDS
                                                ID=rna61088.1.CDS;Parent=rna61088
Chr1 Gnomon
                          43895 43947
                                                ID=rna61088.1.exon; Parent=rna61088
Chr1 Gnomon
                          43895 43947
            exon
Chr1 Gnomon
            start_codon
                         43895 43897
                                                ID=rna61088.1.start_codon;Parent=rna61088
                                           +
                                             1 ID=rna61088.2.CDS;Parent=rna61088
             CDS
                          48839 49007
Chr1 Gnomon
                                          +
                                                ID=rna61088.2.exon; Parent=rna61088
Chr1 Gnomon
                          48839 49007
             exon
                          53889 54000
                                             0 ID=rna61088.3.CDS;Parent=rna61088
             CDS
Chr1 Gnomon
                          53889 54000
                                            . ID=rna61088.3.exon; Parent=rna61088
Chr1 Gnomon
             exon
             CDS
                          55055 55173
                                        . + 2 ID=rna61088.4.CDS;Parent=rna61088
Chr1
     Gnomon
                          55055 55173
                                                 ID=rna61088.4.exon; Parent=rna61088
Chr1
     Gnomon exon
```

GFF3

GFF Header: Pragma begin with "#", comments with "#". Format pragma required for GFF3, highly-recommended for GFF2/GTF.

Pragma/Directives

Pre-declared set of pragma with specific formats/definitions.

Mostly for computers/browsers.

```
##aff-version 3
##spectes http://www.ncbi.nlm.nih.gov/Taxonomy/Browser/wwwtax.cgi?id=436495
##genome-build RexBase Trex1
##sequence-region Chr1 1 217471166
# Note Trex_genome.fasta, complete genome
Chr1 Gnoman gene
                           43895
                                                  ID=gene32251; Name=L0C101732307
ChrFormatiVersion
                           43895
                                                  ID=rna61088; Name=XM_012954515.1; Parent=gene32251
                                                  ID=rna61088.1.CDS;Parent=rna61088
                           43895
  Pragma/Directive
                                  Comments
                           43895
                                                  ID=rna61088.1.exon; Parent=rna61088
 Required for GFF3, don
                           43895
                                                  ID=rna61088.1.start_codon;Parent=rna61088
                                 Free-form text
                           48839
                                                  ID=rna61088.2.CDS; Parent=rna61088
highly-recommended
                                  for humans,
                           48839
                                                  ID=rna61088.2.exon; Parent=rna61088
for GFF2/GTF formats
                           53889
                                                  ID=rna61088.3.CDS; Parent=rna61088
                                  5ignored by ∅
                                                  ID=rna61088.3.exon; Parent=rna61088
                                    parsers.
                                                  ID=rna61088.4.CDS; Parent=rna61088
                                                  ID=rna61088.4.exon; Parent=rna61088
      Gnomon exon
```

GFF3

Feature Attributes

Semi-colon separated

GFF Features: Nine tab-delimited fields required. Null values a "."

```
Feature Strand
                                                                              Key=Value pairs;
      Reference ID
                                                   either "+" or "-",
                                                                          reserved keys begin with
                                Feature Type
Chromosome/scaffold ID
                                                      or "." if no
                                                                               capitals letters;
                             Must be SO term or
    May only contain
                                                    strand applies
                                                                          "Parent" attribute defines
                              accession number
    characters in set:
                                                                         feature hierarchy; must use
 [a=zA-Z0-9::^*$@!+ ?-|1
                                    Score
                                                                              URL-escaping for
                                                        /www.tax.cai?id=436495
                                 floating point
                                                                            forbidden characters
                            complete gend
                                                    ID=gene32251; Name=L0C101732307
                             43895
        Gnomon
                                    78350
                aene
                mRNA
                             43895
                                    78350
                                                    ID=rna61088; Name=XM_012954515.1; Parent=gene32251
  Chr1
        Gnomon
                                                    ID=rna61088.1.CDS; Parent=rna61088
        Gnomon
                CDS
                             43895
                                    43947
  Chr1
                                                    ID=rna61088.1.exon; Parent=rna61088
  Chr1
        Gnomon
                             43895
                                    43947
                exon
                start_codon
                             43895
                                    43897
                                                    ID=rna61088.1.start_codon;Parent=rna61088
        Gnomon
  Chr1
                                                       rna61088.2.CDS; Parent=rna61088
        Gnomon
                                                          61088.2.exon; Parent=rna61088
        Gnomon Source
                                                             88.3.CDS:Parent=rna61088
                                   Start and End
       sually the program or 889
                                                              8.3.exon; Parent=rna61088
                                    55Positions<sup>2</sup>
                                                      Codon Phase arent=rna61088
        organization that
                                                    in either 0, 1; expraparent=rna61088
                                       1-based
    generated the annotations
                                   coordinates on
                                                       Offset to next
                                      "+" strand
                                                      codon position
```

Resources

File manipulation/filtering

| pysam (API) | FASTA/Q, BED,
B/CR/SAM, B/VCF | https://pysam.readthedocs.io/en/latest/api.html#sam-bam-cram-files |
|---------------|----------------------------------|--|
| pyFaidx (API) | FASTA | https://doi.org/10.7287/peerj.preprints.970v1 |
| Seqtk | FASTA/Q | https://github.com/lh3/seqtk |
| Seqkit | FASTA/Q | https://doi.org/10.1371/journal.pone.0163962 |
| seqmagick | Many | https://seqmagick.readthedocs.io |
| bedtools | BAM, BED, GFF, VCF | https://bedtools.readthedocs.io |
| bcftools | VCF/BCF | https://samtools.github.io/bcftools |
| genometools | FASTA/Q, GFF3 | http://genometools.org |
| samtools | FASTA/Q, B/SAM | https://github.com/samtools/samtools |
| vcftools | VCF/BCF | https://vcftools.github.io/man_latest.html |

Resources

QA/QC, Adapter and Quality trimming

| trimmomatic | FASTQ | http://usadellab.org/cms/?page=trimmomatic |
|-------------|--------------|--|
| FastQC | FASTQ, B/SAM | https://www.bioinformatics.babraham.ac.uk/projects/fastqc/ |
| Sickle | FASTA/Q | https://github.com/ucdavis-bioinformatics/sickle |
| Scythe | FASTA/Q | https://github.com/ucdavis-bioinformatics/scythe |
| Sabre | FASTA/Q | https://github.com/najoshi/sabre |
| cutadapt | FASTA/Q | https://cutadapt.readthedocs.io/en/stable/ |

Alignment

| minimap2 | FASTA/Q | https://github.com/lh3/minimap2 |
|-----------|---------|--|
| BWA | FASTA/Q | https://github.com/lh3/bwa |
| bowtie2 | FASTA/Q | http://bowtie-bio.sourceforge.net/bowtie2/index.shtml |
| STAR | FASTQ | https://github.com/alexdobin/STAR |
| GMAP | FASTA/Q | http://research-pub.gene.com/gmap/ |
| exonerate | FASTA | https://www.ebi.ac.uk/about/vertebrate-genomics/software/exonerate |

Resources

Variant calling

| FreeBayes | BAM, VCF | https://github.com/ekg/freebayes |
|-------------|----------------------|--|
| GATK4 | FASTA/Q, B/CRAM, VCF | https://software.broadinstitute.org/gatk/documentation |
| DeepVariant | FASTA/Q | https://github.com/google/deepvariant |
| vg | FASTA/Q | https://github.com/vgteam/vg |

Common file issues

- Non-printable characters
- Non-ASCII encoded characters
- Incorrect formatting (spaces instead of tabs)
- Truncated files

